F. BOSCH
MEDIZINTECHNIK

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FETAL DOPPLER **EN**

FETAL DOPLER **TR**

DOPPLER FŒTAL **FR**

DOPPLER FETAL **ES**

FETALER HERZMONITOR **DE**

Instruction Manual **EN**

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Manual de instrucciones **ES**

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Please read the User Manual carefully to make sure safe and proper use of this Fetal Doppler, Please read and fully understand the Safety Precautions before use.

Precaution Labels Definition

The signal words shown below, left, identify the potential hazard categories. The definition of each category is as follows:

DANGER: This alert identifies hazards that will cause serious personal injury or death.

WARNING: This alert identifies hazards that may cause serious personal injury or death.

CAUTION: This alert identifies hazards that may cause minor personal injury, product damage, or property damage.

NOTE: The label indicates what you should attention.

CONTENT

Section 1 : Introduction 1	
1.1 Contact Information.....	1
1.2 Product Description.....	2
1.3 Operating Principle.....	2
1.4 Indication for Use.....	3
1.5 Contraindications for Use.....	3
1.6 Note for home use.....	3
Section 2 : Safety Guidance.....	3
2.1 Safety Alter Descriptions.....	3
2.2 Symbols.....	5
Section 3 : Getting Started.....	7
3.1 Unpacking an Inspecting.....	7
4.1 Appearance.....	8
4.2 Display Panel.....	9
4.3 Buttons function.....	10
Section 5 : Operation.....	10
5.1 Installation battery.....	10
5.2 Charging the battery.....	11
5.3 Operation instruction.....	11
5.4 Guidelines for manual calculation FHR.....	12
5.5 Modes Setting.....	13
5.6 Detecting.....	13
Section 6 : Maintenance and Services.....	14
6.1 Cleaning.....	14
6.2 Disinfection.....	16
6.3 Maintenance.....	16
6.4 Authorized Repair Service.....	17
Section 7 : Technical Specifications.....	17
Section 8 : Troubleshooting.....	19
Appendix A : EMC information.....	20

Section 1 : Introduction

Become familiar with the controls and how to use the product properly before operating the product.

CAUTION: Federal law restricts this device to be sold by or on the order of a physician or practitioner licensed by state law in which he/she practices to use or order the use of the device.

CAUTION: It should not be used in life supporting or life sustaining applications.

CAUTION: It cannot replace the professional fetal monitor, when the fetal heart rate is abnormal, or can not find the fetal heart, can not feel the fetal movement, pregnant woman's should immediately go to the hospital to seek the doctor's help.

Intended use

The Fetal Doppler Babygeräusch is a hand-held, rechargeable audio Doppler device used for detecting fetal heartbeats.

1.2 Product Description

The Product is a lightweight, portable, detector. It is designed to meet your detecting and hearing needs by providing advanced detecting functions and a full range of sound of the fetal heartbeat.

The product is mainly used to detect the fetal heartbeat rate FHR and the sound of the fetal heartbeat SFH

The growth and development of a fetus can be found out through examination of these indices. It is applicable for department of gynaecology and obstetrics and clinic daily.

In accordance with classification criteria in Annex IX on “Medical Device Directive 93/42/EEC”, the product is class II based on rule 10, “Devices for Direct Diagnosis or Detection on physiological process”.

The Product is powered by 7.2V NI-MH. And Probe of Babygeräusch is can't replace.

1.3 Operating Principle

Fetal Doppler consists of probe (transmitter and receiver) and signal process unit.

Ultrasonic wave is transmitted from one piezoelectric ceramic at the front of the probe to the uterus of the pregnant women. Echo is received by the other piezoelectric ceramic at the front of the probe when ultrasonic wave reaches the fetal heart. Then it is converted into voltage. This Doppler signal is detected and demodulated from the received signal. And the Doppler frequency is consistent with the rhythm of the fetal systole and diastole. Once cardiac valves vibrate and a Doppler frequency excursion is formed. It is transmitted an output signal of cardiac valves vibrating, and it is sent to the loudspeaker for getting a rhythmical sound with the fetal heartbeat. Simultaneously, it is sent to a counter which calculate periods of heartbeat that is fetal heartbeat rate (bpm=beat per minute).

1.4 Indication for Use

The device is ultrasonic fetal heart beat detector, which can detect the Fetal Heart Rate. The built in speaker of the device allows for listen to the fetal heartbeat.

The Product normally is applied to fetus of above 16 weeks growth, difference in pregnant mater.

The normal range of fetal heart rate: 110bpm-160bpm.

Listen to SFH:

Operator can listen to the sound of fetal heartbeat from the headset.

Audio record:

The sound of fetal heartbeats can be recorded by a recorder which is connected with the product. As a safety advisement that can only be connected with a recorder complied with the safety requirements of IEC 60601-1.

1.5 Contraindications for Use

Normally none, as a particular case, please consult your doctor.

1.6 Note for home use

This device cannot replace a professional fetal monitor. If the fetal heart rate is abnormal or cannot be located by using this monitor, the pregnant woman should immediately go to the hospital to seek the doctor's help. If fetal movement is not felt by the pregnant woman, immediately go to the hospital to seek the doctor's help.

Section 2 : Safety Guidance

This unit is internally powered equipment, and the degree of protection against electric shock is type B: 

It means that these person connections will comply with permitted leakage currents, dielectric strengths of IEC/EN60601.

2.1 Safety Alert Descriptions

The following is a list of Product safety alerts that appear in this section and throughout this manual. You must read, understand, and heed these safety alerts before attempting to operate the Product.

DANGER: Fire and Explosion Hazard

Do not operate the Product in the presence of flammable gases to avoid possible explosion or fire hazard.

WARNING: Adjacent and/or Stacked Equipment the product should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Product should be observed to verify normal operation in the configuration in which it will be used.

WARNING: Use only Approved Equipment

Do not use other accessories that are not approved by the manufacturer, such as batteries, gel, cables, or optional equipment. Otherwise it will cause the product function is not normal.

WARNING: Practice the ALARA principle

We recommend that exposure to ultrasound should be kept as low as reasonably achievable principle. This is considered to be good practice and should be observed at all time.

WARNING: Aid healthcare professional tool

The product should not be used in place of normal fetal monitoring. It is a tool to aid the healthcare professional.

CAUTION: Temperature/Humidity/Pressure Extremes Exposing the Product to extreme environmental conditions outside of its operating parameters may compromise the ability of the Product to function properly.

CAUTION: Battery Disposal

Recycle or dispose of the battery in accordance with all federal, state and local laws. To avoid fire and explosion hazard, do not burn or incinerate the battery.

CAUTION: Possible Radio Frequency (RF) Susceptibility

RF susceptibility from cellular phones, CB radios and FM 2-way radio may cause interference with the Product. Do not operate wireless radiotelephones in the vicinity of the Product – turn power OFF to the radiotelephone and other like equipment near the Product.

**CAUTION: Systems Statement**

Equipment connected to the Product must be certified to the respective IEC Standards (i.e. IEC 60950 for data processing equipment and IEC 60601-1 for medical equipment).

**CAUTION: Case Cleaning Solutions**

When disinfecting the case, use a non-oxidizing disinfectant, such as ammonium salts or a glutaraldehyde based cleaning solution, to avoid damage to the metal connectors.

**CAUTION: Environment of use**

The Product is designed for indoor use. Operator must confirm that the environment of use meets the required operating environmental specifications before using.

**CAUTION: Cold Environments**

If the Product is stored in an environment with a temperature below the operating temperature, the unit should be allowed to warm up to the needed operating temperature before using.

**CAUTION: Product recycling or disposal**



The device could be sent back to the manufacturer for recycling or disposal comply with local laws when equipment and accessories at the end of their useful lives.

2.2 Symbols

The following symbols may appear in this manual, on the Product, or on its accessories. Some of the symbols represent standards and compliances associated with the Product and its use.

Consult instructions for use of the Product and/or its accessories.

	Warning Information.
	Authorized Representative in the European Community
	CE Mark: The Product system conforms to essential requirements of the Medical Device Directive 93/42/EEC.
	Date of manufacture.
	Manufacturer Information.
	Storage and transport temperature limit
	Humidity limitation
	Atmospheric Pressure limitation
	CLASS II equipment
	Type B applied part.
<p data-bbox="184 1243 244 1269">IPX4</p>	Ingress Protection.

	Specifies serial number of the Product.
	Production batch.
	It indicates that the equipment should be sent to the special agencies according to local regulation for separate collection after its useful life.
	The environmental protection use period is 5 years.

Section 3 : Getting Started

This section presents information on unpacking and setting up the Product

3.1 Unpacking an Inspecting

Every attempt is made to ensure your order is accurate and complete. However, to be sure that your order is correct, verify the contents of the box against your packing slip.

The Product is designed for simplicity of operation and set-up and requires minimal assembly.

The accessories for normal use supplied or approved by Trimpeks can be used with Babygeräusch.

The following item is packing list:


No.	Item	Quantity	Packed or Unpacked
1	Packing List	1	√
2	Instruction Manual	1	√
3	mainframe(including Probe)	1	√
4	adapter	1	√

Carefully inspect each item as it is unpacked for any signs of damage which may have occurred during shipment.

- Check the components according to the packing list.
- Check for any damage or defects. Do not attempt to setup the Product if anything is damaged or defective. Contact Trimpeks.

Section 4 : Appearance and Structure


4.1 Appearance

- ① LCD display
- ② Volume down button
- ③ Power on/off button
- ④ Loudspeaker
- ⑤ Andio Output Socket
- ⑥ Charging socket

- ⑦ Ultrasonic probe
- ⑧ Volume up button
- ⑨ Mode button
- ⑩ Battery cover
- ⑪ Label

4.2 Display Panel

FHR refresh frequency:

When monitoring, the heart shape symbol will appear and disappear. When it appeared or disappeared one time, the fetal heart beats one time.

FHR Display :

Here displays the fetal heart rate value, whose unit is bpm(beat per minute)

Warning:

When open the lower and higher restrict warn sound, and the FHR is in the restrict scope, the symbol will twinkling. And the warn sound will sound.

Average FHR mode:

When the symbol is appeared, the working mode is Average FHR mode.

Transient FHR mode:

When the symbol is appeared, the working mode is Transient FHR mode.

Battery Indicator:

Display the quantity of electricity in batteries now. When the four black small rectangles are all appeared, it responses batteries are fulfill. As decreases of the quantity of electricity in batteries, the rectangles will disappear from left to right gradually.

Volume Level :

Display the Volume Level.

4.3 Buttons function


- ① **Power on/off button** “” **Function:** power on or off the Doppler
- ② **Volume up button** “+”
Function: increase the volume.(more detail refer to 5.3 Modes Setting)
- ③ **Volume down button** “-”
Function: decrease the volume. (more detail refer to 5.3 Modes Setting)
- ④ **MODE button** “MODE”
Function: switch working mode. (more detail refer to 5.3 Modes Setting)

Section 5 : Operation

This section provides the description for operation.

5.1 Installation battery

Power of the Doppler is a 7.2V NI-MH which is supplied from Trimpeks. The battery was installed by the manufacturer.

WARNING: Don't throw batteries in fire because this may cause it explode.

CAUTION: The battery must be proper disposed according to local regulation after their use

CAUTION: Don't replace rechargeable batteries with ordinary dry batteries.

5.2 Charging the battery

When the battery is lower, Battery Indicator will flicker. You have to shut down the Doppler to charge the battery.

Insert the DC plug of the special charger supplied by Trimpeks in to the charging socket on the Doppler, then connecting the AC plug of charger to the AC power socket.

This process will spend 10 hours to charge the battery. When charging, the Battery Indicator will display and round with a special regular. When the battery is fully, the Battery Indicator is stillness .

WARNING: Charge the battery with the special charger supplied by Trimpeks.

CAUTION: You must shutdown the doppler before charging the battery.

CAUTION: Don't use the doppler for monitoring until the adapter disconnected with the doppler.

5.3 Operation instruction

Operational guidelines

1. Turn on

Press and hold the Power on button for 2 seconds, the device will turn on and self check, when the FHR display "---" , the self-check is finished and get into use.

2. Detect FHR

Put the gel on the surface of the ultrasonic probe, then press the probe on the detect position softly, and ensure it has a good touch surface. When the probe is on the right position, you can hear the fetal heart beat. Wait until the sounds is clear and consecutive seconds later, the LCD screen will display the right FHR.

3. Volume Control

Under the FHR display state, you can control the volume by volume up button and volume down button. Volume level has 9 level , "0" is silent , "8" is loudest .

4. Turn off

After detecting, press and hold the power off button for 2 seconds, the device will turn off .The device will power off automatically without signal one minute.

5.4 Guidelines for manual calculation FHR

1. After turn on the Doppler, press " MODE" button and enter the interface for manual calculation FHR.

2. Detect the FHR

Put the gel on the surface of the ultrasonic probe, then press the probe on the detect position softly, and ensure it has a good touch surface. When the probe is on the right position, you can hear the fetal heart beat. Press down "+" or "-" button, calculate 15 heart beats, then press once more "+" or "-", the screen will display fetal heart rate. If you continue to detect, you can repeat.

3. Switch to Automatic calculation

Press "MODE" button four times

4. Turn off

After detecting, press power off button for 2 seconds, the device will turn off. The device will power off automatically without signal one minute.

5.5 Modes Setting

Function Menu:

Turn on the Doppler. Press the "MODE" button for function set. Operation as follows:

Press first time:

Enter the function of "Manual Calculate FHR", the Screen flickers "000". Press once "+" button or "-" button and calculate 15 heart beats, then press once more "+" or "-" button, the screen will display Fetal Heart rate.

Press second time:

Enter the interface of FHR display mode options, the screen flickers "A" or "T". "A" display the Average Fetal Heart Rate, "T" display Transient Fetal Heart Rate. Default factory settings is: "T". press "+" or "-" button, select the "A" Mode or "T" Mode.

Press third time:

Enter the interface of backlight switch settings, and  is flicking, press "+" or "-" button, select the switch backlight, Default factory settings is: "ON"

Press fourth time:

Enter the alarm setting interface, and  (Upper left corner) is flicking. Press "+" or "-" button select to turn on or off the alarm. Default factory setting is: "OFF".


Notice: This interface is only to set the alarm on or off, so you can't detect the fetal heart rate.

5.6 Detecting

Locate the position of the fetus by hand touching, firstly to find out the best direction to the fetal heart. Place the faceplate of probe at the best position for detecting fetal heart beat. Adjust the probe to obtain an optimum audio signal ideally by angling the probe around. Generally, the site of heart of fetus is 1/3 below of navel line at its earlier stage, it then moves upward with increasing of gestational period, and the site of heart of fetus will be a

little deviation to left or right with different fetus. please make sure that the surface of probe should be contacted fully with the skin. After the sound become clear and stays stable for a few seconds, the FHR value will appear on the LCD in real time. If no coupling gel, water can be used.

Note: abnormal value of FHR may appear during the searching of fetal heart.

Section 6 : Maintenance and Services

6.1 Cleaning

Listed below are recommendations for cleaning the Product and its accessories.

Recommended cleaning product:

The following cleaning products may be used to clean the exterior surfaces of the Product as well as the batteries.

- Isopropyl alcohol (70% solution in water)
- Mild soap and water
- Sodium hypochlorite (chlorine bleach) (3% solution in water)
- Quaternary ammonium compounds (such as Lysol) (10% solution in water)

- Do not use abrasive cleaners or strong solvents such as acetone or acetone-based cleaners.
- Do not mix disinfecting solutions (such as bleach and ammonia) as hazardous gases may result.
- Do not clean electrical contacts or connectors with bleach.

Cleaning instructions:

1. Before cleaning the Product, turn off the device and disconnect the power cord.
2. Before cleaning, remove all adherent soil (tissue, fluids, etc.) and wipe thoroughly with a cloth dampened with water before applying the cleaning solution.
3. When cleaning, do not immerse. Keep the exterior surface of the device clean and free of dust and dirt ,clean exterior surface of the unit with a dry, soft cloth .if necessary, clean it with a soft cloth soaked in a solution of soap and wipe dry with a clean cloth immediately.
Wipe the transducer body with soft cloth to remove any remaining coupling gel .Clean with soap only
4. Wring any excess moisture from the cloth before cleaning.
5. Avoid pouring fluids on the device, and do not allow fluids to penetrate the exterior surfaces of the device.
6. To prevent scratching the display, the use of a soft cloth is recommended.

CAUTION: To prevent damage to equipment, do not clean any part of the Product or Accessories with phenolic compounds. Do not use abrasive or flammable cleaning agents. Do not steam, autoclave, or gas-sterilize the product or accessories.

CAUTION: Cleaning liquids: Do not submerge the device in liquids or pour cleaning liquids over, into or onto the device.

- * Don't use strong solvent, for example. Acetone.
- * Never use an abrasive such as steel wool or metal polish.
- * Do not allow any liquid to enter the product, and do not immerse any parts of the device into and liquids.
- * Avoid pouring liquids on the device while cleaning.
- * Don't remain any cleaning solution on the surface of the device.
- * Wipe the surface of sensor of transducer with 70% ethanol or alcohol, self-air dry or clean with a clean, dry cloth.

6.2 Disinfection

Cleaning the unit surface and the transducer as the above mentioned, then wipe the surface of transducer with 70% ethanol or alcohol, clean the transducer surface with a dry, soft cloth.

WARNING: Don't use low temperature steam sterilization or other way to sterilize

WARNING: Don't use high temperature sterilizing process

6.3 Maintenance

WARNING: Failure on the part of all responsible individuals, hospitals or institutions, employing the use of Product, to implement the recommended maintenance schedule may cause equipment failure and possible health hazards. The manufacturer does not, in any manner, assume the responsibility for performing the recommended maintenance schedule. The sole responsibility rests with the individuals, hospitals, or institutions utilizing the product.

1. The transducer acoustic surface is frangible and must be handle with care .Gel must be wiped off from the transducer after use. These precautions will prolong the life of the unit.
The user must check that the equipment does not have visible evidence of damage that may affect patient safety or product's capability before use .The recommended inspection interval is once

per month or less. If damage is evident, replacement is recommended before use.

2. To ensure the product is always functional when required, the following maintenance shall be performed:

- Visual Inspection
- Cleaning the product and its accessories
- Check the battery power
- Testing product Performance

Correction: manually calculate the FHR with hearing fetal heartbeat sound for qualification.

Recommended maintenance and care:

- It is important that the Product is stored at the operating temperature range if it is expected to be used. Optimal battery life will be obtained if stored and operated at room temperature. See Section 7 for Temperature Specifications.
- The Product requires no calibration.

6.4 Authorized Repair Service

The Product has no user-serviceable internal components. Try to resolve any maintenance issues with the Product by using the Troubleshooting Table presented in this chapter. If you are unable to resolve the problem, contact Trimpeks.

NOTE: The warranty will be void upon unauthorized disassembly or service of the Product.

Section 7 : Technical Specifications

This section presents the specifications and safety standards of the product.

NOTE: The following specifications are subject to change and are only noted as a point of reference.

Product name: Fetal Doppler

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Safety: Complies with IEC 60601-1, EN60601-1-2, IEC 61266

Physical Dimensions:

135mm W x 83mm D x 33 mm H, Wt: 0.24kg (not including battery)

Classification:

Anti-electric Shock Type: Internally powered equipment

Protection against electric shock: CLASS II equipment 

Anti-electric Shock Degree: Type B equipment 

Degree of Protection against Harmful Ingress of Water

Main Unit : Non-protected

Probe: IPX4 Water Ingress Protection Code

Degree of Safety IN Presence of Flammable Gases: Equipment not suitable for use in presence of flammable gases

Working Mode: Continuous operation

Technical parameters:

Ultrasound:

Ultrasonic emitting frequency: 3MHz

Overall sensitivity at the distances 200mm from the face of the probe (Doppler frequency: 500 ± 50 Hz, Target velocity: 10cm/s~40cm/s): 92.9dB

Spatial-peak temporal-peak acoustic pressure: 20.8KPa Output power: 8.6mW

Effective area of the ultrasonic transducer active element: 3.07cm²

The acoustic coupling medium for normal use: ph: 5.5~8, Acoustic impedance: $\leq 1.7 \cdot 10^5$ g/cm²·s

Working mode: Continuous wave Doppler Audio output:

Audio output power: <0.5 W

Audio out Socket: $\Phi 3.5$ mm

Product service life: 5 years

Recommended battery type:

7.2V NI-MH, which is supplied from Trimpeks.

Adapter:The input voltage see the label

Charge limit: DC 15V 1.5A(max)

Stand-by Time: >8hours

Performance:

Sensitivity: 16 weeks gestation

FHR Measuring Range: 50-210bpm (beat per minute Resolution:

1bpm

Accuracy: ± 2 bpm

Environmental Requirements:

Operating Conditions:

Temperature: 5° C to 40° C (41°F-104°F)

Humidity: <80% RH, non-condensing

Atmospheric pressure: 86kPa to 106kPa

Storage and Shipping Conditions: Temperature: -

20° C to 55° C (-4°F-131°F)

Humidity: 10% - 93% RH, non-condensing

Atmospheric pressure: 50kPa to 106kPa

Section 8 : Troubleshooting

What appearance to be a malfunction may not always serious, if your device should not perform as expected, consult the table below to see if the problems can be corrected before seeking help from your dealer or service representative.

Symptom	Cause	Remedy
Audio scream	<ul style="list-style-type: none"> ● Audio volume too high ● Too much gel on the probe surface ● Battery is exhausted 	<ul style="list-style-type: none"> ● Turn down the volume ● Use less gel ● Change or charging the battery
Weak sound	<ul style="list-style-type: none"> ● Audio volume too low ● Insufficient gel 	<ul style="list-style-type: none"> ● Turn up the volume ● Add gel

output	<ul style="list-style-type: none"> ● Battery is exhausted 	<ul style="list-style-type: none"> ● Charging the battery
Low sensitivity	<ul style="list-style-type: none"> ● Incorrect probe position ● Insufficient gel 	<ul style="list-style-type: none"> ● Locate the correct position ● Add gel

Appendix A : EMC information

CAUTION: Fetal Doppler needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided for in the ACCOMPANYING DOCUMENTS.

CAUTION: Portable and mobile RF communications equipment can affect **Fetal Doppler**.

CAUTION: the **Fetal Doppler** should not be used adjacent to or stacked with other equipment.

Guidance and manufacturer's declaration-electromagnetic emissions		
The models Babygeräusch are intended for use in the electromagnetic environment specified below. THE customer or the user of the models Babygeräusch should assure that it is used in such an environment.		
Emissions	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The models Babygeräusch use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The models Babygeräusch is suitable for use in all establishments, including domestic establishments and those directly connected to the public
Harmonic emissions IEC 61000-3-2	Class A	
Voltage	Not applicable	

fluctuations/flicker emissions IEC 61000-3-3		low-voltage power supply network that supplies buildings used for domestic purposes.	
Guidance & Declaration—electromagnetic immunity			
The models Babygeräusch are intended for use in the electromagnetic environment specified below. THE customer or the user of the models Babygeräusch should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for Input/output lines	±2kV for power supply lines	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC 61000-4-5	±1kV line to line ±2kV line to earth	±2kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95 % dip in U_T)for 0.5 cycle 40% U_T (60% dip in U_T)for 5cycles 70% U_T (30% dip in	<5% U_T (>95 % dip in U_T)for 0.5 cycle 40% U_T (60% dip in U_T)for 5cycles 70% U_T (30% dip in U_T)for 25	Mains power quality should be that of a typical commercial or hospital environment. If the user of the models Babygeräusch requires

	U _T)for 25 cycles <5% U _T (>95% dip in U _T) for 5sec	cycles <5% U _T (>95% dip in U _T) for 5sec	continued operation during power mains interruptions, it is recommended that the models Babygeräusch be powered from an uninterruptible power supply or a battery.
Power frequency(50/60 Hz)magnetic field IEC 61000-4-8	3 A/m	Not applicable	Not applicable
NOTE U _T is the a.c. mains voltage prior to application of the test level.			

The models Babygeräusch is intended for use in the electromagnetic environment specified below. The customer or the user of the Fetal Heart Rate Doppler Detector should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 Vrms 80 MHz to 2.5 GHz	3V 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the models Babygeräusch including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance

			<p>3V</p> <p>$d=1.2 \times P^{1/2}$ 80MHz to 800 MHz</p> <p>$d=1.2 \times P^{1/2}$ 800MHz to 2.5 GHz Where</p> <p>P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the models Babygeräusch is used exceeds the applicable RF compliance level above, the models Babygeräusch should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the models Babygeräusch.
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the models Babygeräusch.

The models Babygeräusch is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the models Babygeräusch can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the models Babygeräusch as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.2 \times P^{1/2}$	80 MHz to 800 MHz $d=1.2 \times P^{1/2}$	800 MHz to 2,5 GHz $d=2.3 \times P^{1/2}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the

transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Overall sensitivity(3MHz Probe)

Temperature: 24.9°C, Relative Humidity : 68 %

Distance for the probe face	50mm
-----------------------------	------

Nominal acoustic working frequency	3MHz
------------------------------------	------

Doppler frequency	296Hz
-------------------	-------

Target Velocity	12cm/s
-----------------	--------

A(d): Target reflection loss	44.5dB
------------------------------	--------

B: two-way attenuation over the acoustic pathway	52.5dB
--	--------

C: Signal-to-noise Ratio	Vs: 298mV Vn: 119mB $C=20\log(Vs/Vn)=8.0dB$
--------------------------	---

S: Overall Sensitivity	$S=A+B+C=105db$
------------------------	-----------------

Supplementary information: --

Distance for the probe face	75mm
-----------------------------	------

Nominal acoustic working frequency	3MHz
------------------------------------	------

Doppler frequency	296Hz
-------------------	-------

Target Velocity	12cm/s
-----------------	--------

A(d): Target reflection loss	44.5dB
------------------------------	--------

B: two-way attenuation over the acoustic pathway	46.2dB
--	--------

C: Signal-to-noise Ratio	Vs: 293mV Vn: 117mB $C=20\log(Vs/Vn)=8.0dB$
--------------------------	---

S: Overall Sensitivity	$S=A+B+C=98.7db$
------------------------	------------------

Supplementary information: --

Distance for the probe face	100mm
Nominal acoustic working frequency	3MHz
Doppler frequency	296Hz
Target Velocity	12cm/s
A(d): Target reflection loss	44.5dB
B: two-way attenuation over the acoustic pathway	45.1dB
C: Signal-to-noise Ratio	Vs: 291mV Vn: 116mB $C=20\log(Vs/Vn)=8.0\text{dB}$
S: Overall Sensitivity	$S=A+B+C=97.6\text{db}$
Supplementary information: --	

Distance for the probe face	200mm
Nominal acoustic working frequency	3MHz
Doppler frequency	296Hz
Target Velocity	12cm/s
A(d): Target reflection loss	44.5dB
B: two-way attenuation over the acoustic pathway	40.2dB
C: Signal-to-noise Ratio	Vs: 289mV Vn: 113mB $C=20\log(Vs/Vn)=8.2\text{dB}$
S: Overall Sensitivity	$S=A+B+C=92.9\text{db}$
Supplementary information: --	

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**Shenzhen Jumper Medical Equipment Co.,Ltd**
MedPath GmbH