

GENTLE

WRIST-TYPE FULLY AUTOMATIC DIGITAL BLOOD PRESSURE MONITOR

- 1 Average of last 3 measurements
- 2 Accurate measurement
- 3 Irregular heartbeat detection
- 4 Quick Result
- 5 Quiet inflation method
- 6 Automatic power-off
- 7 WHO Classification indicator



SPECIFICATIONS

Product Description	Wrist-type Fully Automatic Blood Pressure Monitor	
Model	F.BOSCH GENTLE	
Display	LCD Digital Display Size:45mm x 30.5mm(1.77" x 1.20")	
Measurement Method	Oscillometric Method	
Measurement Range	Systolic Pressure	60mmHg~260mmHg
	Diastolic Pressure	40mmHg~200mmHg
	Pressure	0mmHg~299mmHg
	Pressure	±3mmHg
	Pulse	30 ~ 180 Beats/Minute
	Pulse	±5%
Pressurization	Automatic Pressurization	
Memory	2 x 150 Memories in Two Groups with Date and Time	
Function	Irregular Heartbeat Detection	
	WHO Classification Indicator	
	Last 3 Results Average	
	Low Battery Detection	
	Automatic Power-Off	

Power Source	2 Alkaline Batteries Size AAA	
Battery Life	Approximately 2 months at 3 tests per day	
Unit Weight	Approx. 91g (2.21 oz) (Excluding Battery)	
Unit Dimensions	Approx. 82.3mm×63.8mm×27.5mm(L x W x H) (3.24" x 2.51" x 1.08")	
Cuff Circumference	Fits wrist circumference 13.5-21.5 cm(5.3"-8.5")	
Operating Environment	Temperature	10°C ~ 40°C (50°F~104°F)
	Humidity	15%~93%RH
	Pressure	800hPa~1060hPa
Storage Environment	Temperature	-25°C~55°C (-13°F~131°F)
	Humidity	≤93% RH
Ingress Protection Rating	IP 22	
Classification	Internal Powered Equipment Type BF 	