MOBINEB MESH NEBULIZER



- Rechargeable battery operated
- 2 Replaceable medication cup
- 3 Two levels nebulization rate
- 4 Particle size <5 µm
- 5 Big medication cup

- 6 Low noise
- Light-weight
- 8 Dual power supply
- 9 Accessories:
 - · Medication cup · Child mask
 - · Adult mask
- · USB cable

FEATURES	
Power supply:	DC 5V or (3.7V rechargeable lithium battery)
Input power:	<15 VA
Nebulization rate in level I:	≥ 0.25mL/min
Nebulization rate in level II:	≥ 0.15mL/min
Noise:	≤50 dB
Weight:	0.1 kg
Equivalent volume particle diameter distribution:	The volume distribution ratio of small atomized particles (diameter $\!<\!5~\mu m$) is no less than 60 $\%$
Type of protection against electric shock:	Class II equipment, internally powered equipment
Degree of protection against electric shock:	Type BF applied part
Degree of protection against ingress of liquid:	IP22

Note: Please purchase medical power adapter from qualified manufacturer



OPERATIONAL ENVIRONMENT

(input: AC 100-240V, 50Hz-60Hz, output: DC 5V, 1.0A)

Temperature: 5 °C \sim 40 °C Humidity: 15 % \sim 90 %

Atmospheric pressure: 700 hPa \sim 1060 hPa

TRANSPORT AND STORAGE

Environment of transport and storage:

Temperature: -40 °C \sim +55 °C Relative humidity: 5 % \sim 96 %

Atmospheric pressure: 500 hPa ~ 1060 hPa

